

Product update: HEINE® LoupeLight 2

With HEINE HR or HRP Binocular Loupes for a perfect view.

The combination of light and loupe delivers the ultimate in performance as high quality illumination with LED^{HQ} maximises the image quality of the loupe.

For a perfect view and a precise diagnosis, the LoupeLight 2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight 2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight 2 (product description see overleaf).

NEW!

- ✓ Now 25% brighter due to improved LED technology
- ✓ Replaceable connecting cord



Replaceable connecting cord with high-quality splashproof threaded connector for heavy duty use.

Now even 25% brighter

LED^{HQ}

LED NOW IN HEINE QUALITY.

HR 2.5x High Resolution Binocular Loupes. Achromatic optics for crisp, bright images.

Product update: HEINE® MicroLight 2

Ultralight, universal headlight.

For those who need a compact, lightweight illumination solution without magnification, we offer the MicroLight 2. Featuring LED^{HQ} illumination for homogeneous illumination with true and natural colour rendering for every examination situation. Two mounting options offer individuality and easy handling. Illumination without compromises, in a lightweight design.



HEINE MicroLight 2 on S-FRAME



HEINE MicroLight 2 on Lightweight Headband

See overleaf for ordering information.



HEINE® mPack mini – The compact and light power source for LoupeLight 2 and MicroLight 2.

The new, high performance power source HEINE mPack mini for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of just 95g and dimensions of 44 x 23 x 101 mm, this mobile power source is now approximately approx. 70% smaller and lighter than its predecessor.

The pouch, which is included in the delivery, can be attached to the belt. The mPack mini can also be charged whilst in the pouch. Thanks to the Li-ion battery technology the power source has an operating time of typ. 9 hours, in continuous operation at maximum power typ. 4 hours. The mPack mini is flexible in charging due to its USB connector. Charge it via laptop, desktop or USB charger.



 **USB charging option**

NEW!

- ✓ 70% smaller and lighter
- ✓ Operating time typ. 9 hours
- ✓ Charging time typ. 4 hours



Size comparison of HEINE mPack LL/ mPack mini

LoupeLight 2 Sets with HR Binocular Loupes

HR 2.5x optics with i-View on S-FRAME, LoupeLight 2, detachable yellow filter, mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid, cleaning cloth, carrying case

340mm working distance	C-008.32.450
420mm working distance	C-008.32.451
520mm working distance	C-008.32.452

MicroLight 2 Sets

MicroLight 2 on S-FRAME with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, retaining cord and cleaning fluid in case


J-008.31.276

MicroLight 2 on Lightweight Headband with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, case

J-008.31.277

Single Devices

LoupeLight 2 for HR/HRP Loupes with i-View loupe mount	C-008.32.226
mPack mini with Li-ion rechargeable battery and E4-USB plug-in transformer	X-007.99.650
mPack mini with Li-ion rechargeable battery, without E4-USB plug-in transformer	X-007.99.649
Connecting cord SC 1 (1.6m) for LoupeLight 2 or MicroLight 2/ mPack mini	C-000.32.240

For further ordering information please contact your HEINE distributor 


 7248 S. Tucson Way
 Centennial, CO 80112
 USA Headquarters
 Heine USA


 AUTHORIZED DEALER

