



# ELITE WIRELESS LED HEADLIGHT

FINEST QUALITY IN THE INDUSTRY / AGGRESSIVE PRICING



Wrap around rear cranial support on headband for added comfort. Removable if needed



Removable / sterilizable lever to safely adjust beam direction

Rheostat for spot variability 20-110 mm at 16"/40 cm working distance

Ideal color temperature at 4,500° K

 Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

**Sontec Instruments Elite Wireless LED Headlight offers the ability to be untethered from a light source or remote battery pack, while providing bright, white light with a homogenous spot. The lightweight RCS headband provides comfort and stability where a wireless headlight is preferred**

## Features

- 10w specialized LED offers the BRIGHTEST illumination available
- Light Intensity: 70,000 Lux @ 16"/40cm
- Module with variable spot size: 20mm-110mm @ 16"/40cm
- Lightweight battery and clip on holster provide over 3 hours of continuous light at maximum intensity
- Variety of headband options available for ultimate comfort
- Audible low battery warning

## Specifications

- 10w LED Producing Bright, White Light
- Color temperature: 4,500°K
- Typical LED life: >50,000 hours
- Lightweight Headband with Rear Cranial Support for added comfort
- Variable Module Spot 20mm-110mm @ 16" / 40cm
- Variable Intensity Control: 0 – 70,000 Lux at 16"/ 40cm working distance
- Battery Life – 3+ and 6+ hour batteries available
- Battery Charger Accepts 100-240V AC, 50/60 Hz input for worldwide operation
- Battery Recharge Time – 3.5 hrs
- Warranty Three (3) Years on LED Headlight and One (1) Year on Battery

## Classification

- European Class I, Active device per Annex IX, rule 1
- FDA Class I, Device Listing 125713
- Complies with Essential Requirements Matrix of MDD 93/42/EEC, amended by 2007/47/CE

## System Includes

- LED Headlight with RCS Headband
- Battery Pack with 2 Rechargeable Batteries
- Single or Dual Bay Charger
- Padded Carrying Bag
- Joystick for positioning Module

