



 Sontec Instruments®

KERRISON RONGEUR


Kerrison Universal Rotating 360° Rongeur

Quick trigger release / shaft locks in 8 positions through a 360° rotation / 8" (20cm) shaft

- 404-606 handle
- 404-605 full set includes: handle all 40° sizes 8", all 90° sizes 8"
- 404-6071 1 mm bite 90°
- 404-6072 2 mm bite 90°
- 404-6073 3 mm bite 90°
- 404-6074 4 mm bite 90°
- 404-6075 5 mm bite 90°
- 404-6081 1 mm bite 40°
- 404-6082 2 mm bite 40°
- 404-6082B 2 mm bite 40° up curved
- 404-6082C 2 mm bite 40° down curved
- 404-60825E 2.5 mm bit 40°
- 404-6083 3 mm bite 40°
- 404-6083B 3 mm bite 40° up curved
- 404-6083E 3 mm bite 40° endoscopic
- 404-6084 4 mm bite 40°
- 404-6085 5 mm bite 40°

10" (25.5cm) shaft

- 404-606 handle
- 404-607 full set includes: handle all 40° sizes 10", all 90° sizes 10"
- 404-6076 1 mm bite 90°
- 404-6077 2 mm bite 90°
- 404-6078 3 mm bite 90°
- 404-6079 4 mm bite 90°
- 404-6080 5 mm bite 90°
- 404-6086 1 mm bite 40°
- 404-6087 2 mm bite 40°
- 404-6088 3 mm bite 40°
- 404-6089 4 mm bite 40°
- 404-6090 5 mm bite 40°



Handle sold separately from shaft

Kerrison Universal Rotating 360° Rongeur

404-606B handle *special order

Bayonet shaft / black ceramic / 12 1/2" (31.8 cm)


404-6061B 1 mm bite 40°

404-6062B 2 mm bite 40°

404-6063B 3 mm bite 40°

Endoscopic / black ceramic / 16" (40.6 cm)

404-6063EB 3 mm bite 40°



Sontec Beatty Foraminotomy Rongeur

7 3/4" (19.5 cm) shaft / curved up

- 404-461 2 mm bite
- 404-462 3 mm bite
- 404-463 4 mm bite

7 3/4" (19.5 cm) shaft / curved down

- 404-450 2 mm bite
- 404-451 3 mm bite
- 404-452 4 mm bite


* Also available with rounded footplate



slightly curved up



curved up



curved down




Beatty Foraminotomy Rongeur

Black ceramic coating / 8" (20.3 cm) shaft /
slightly curved up

- 404-461B 2 mm bite
- 404-462B 3 mm bite
- 404-463B 4 mm bite

* Also available with rounded footplate



Kerrison 40° Rongeur


Black ceramic coated / thin
footplate punch / 10 mm opening

		Bite Width	Length of Shaft		
			6" / 15 cm	7" / 18 cm	8" / 21 cm
40° Up-Angle		1 mm	404-490B	404-539B	404-400B
		2 mm	404-491B	404-540B	404-402B
		3 mm	404-492B	404-541B	404-404B
		4 mm	404-493B	404-542B	404-406B
		5 mm	404-494B	404-543B	404-408B
		6 mm	404-495B	404-544B	404-410B



Kerrison 90° Rongeur


Black ceramic coated / 10 mm opening

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
Up-Bite 	1 mm	404-470B	404-519B	404-412B
	2 mm	404-471B	404-520B	404-414B
	3 mm	404-472B	404-521B	404-416B
	4 mm	404-473B	404-522B	404-418B
	5 mm	404-474B	404-523B	404-420B
	6 mm	404-475B	404-524B	404-422B



Kerrison 90° Rongeur


Thin footplate-punch / ceramic coating


	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
 <p>90° Down-Bite</p>	1 mm	404-480B	404-530B	404-424B
	2 mm	404-481B	404-525B	404-426B
	3 mm	404-482B	404-526B	404-428B
	4 mm	404-483B	404-527B	404-430B
	5 mm	404-484B	404-528B	404-432B
	6 mm	404-485B	404-529B	404-434B



Kerrison 40° Laminectomy/Lumbar Rongeur


Black ceramic coating / 8" (20 cm) shaft / wide mouth 15 mm bite
/ with & without ejection pin

	Bite Width	With Ejector Pin	Without Ejector Pin
40° Up-Angle 	1 mm		404-337B
	2 mm	404-332B	404-338B
	3 mm	404-333B	404-339B
	4 mm	404-334B	404-340B
	5 mm	404-335B	404-341B



Kerrison 40° Rongeur


Thin footplate-punch / 10 mm opening

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
 40° Up-Angle	1 mm	404-490	404-539	404-400
	2 mm	404-491	404-540	404-402
	3 mm	404-492	404-541	404-404
	4 mm	404-493	404-542	404-406
	5 mm	404-494	404-543	404-408
	6 mm	404-495	404-544	404-410



Kerrison 90° Rongeur

Thin footplate punch / 10 mm opening / up

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
	1 mm	404-470	404-519	404-412
	2 mm	404-471	404-520	404-414
	3 mm	404-472	404-521	404-416
	4 mm	404-473	404-522	404-418
	5 mm	404-474	404-523	404-420
	6 mm	404-475	404-524	404-422



Sontec Kerrison Rongeur 90° Down Bite Rongeur







Thin footplate / 10mm opening down

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
 90° Down Bite	1 mm	404-480	404-530	404-424
	2 mm	404-481	404-525	404-426
	3 mm	404-482	404-526	404-428
	4 mm	404-483	404-527	404-430
	5 mm	404-484	404-528	404-432
	6 mm	404-485	404-529	404-434



Sontec Kerrison 40° Rongeur


Thin footplate / 10 mm opening / down

		Bite Width	Length of Shaft	
			7" / 18 cm	8" / 21 cm
40° Forward Down-Bite		1 mm	404-5301	404-553
		2 mm	404-5302	404-554
		3 mm	404-5303	404-556
		4 mm	404-5304	404-557
		5 mm	404-5305	404-558
		6 mm	404-5306	404-559



Kerrison 40° Rongeur







Black and gold ceramic coating
/ thin footplate-punch / up / 8"
(20.3 cm)

	Bite Width	Metal
		Black/Gold
	1 mm	404-400BG
	2 mm	404-402BG
	3 mm	404-403BG
	4 mm	404-404BG
	5 mm	404-405BG



Kerrison Bayonet 40° Rongeur


Thin footplate / 10 mm opening / 8" (20 cm) shaft / with & without ejection pin

		Bite Width	With Pin	Without Pin
40° Up-Angle		1 mm	NA	404-401
		2 mm	404-4052	404-409
		3 mm	404-405	404-4033
		4 mm	404-4054	404-407
		5 mm	404-4055	404-403
				



Kerrison Bayonet 40° Rongeur

Black ceramic coating / thin footplate / 10 mm opening / 8" (20 cm) shaft / with & without ejection pin


	Bite Width	With Pin	Without Pin
	1 mm	NA	404-401B
	2 mm	404-4052B	404-409B
	3 mm	404-405B	404-4033B
	4 mm	404-4054B	404-407B
	5 mm	404-4055B	404-403B



Kerrison Bayonet Foraminotomy Rongeur

2 mm opening / 10" (25.4 cm) shaft








404-468



Kerrison Convertible 40° Rongeur

8" (20.3 cm) shaft / 10 mm opening / without
ejection pin / ejector pin for 3 mm sizes and
longer only

- * Other tip sizes & lengths available upon request.
- * Swing system for easy cleaning
- * Sizes available 4-6 mm with pin

		Bite Width		
40° Up-Angle		1 mm	404-337C	without pin
		2 mm	404-338C	without pin
		3 mm	404-339C	with pin
		4 mm	404-340C	with pin
		5 mm	404-341C	with pin
		6 mm	404-342C	with pin
				



Kerrison 40° Rongeur Titanium


7" (18 cm)

2603-541 3 mm



Sontec Kerrison 40° Rongeur

10 mm opening / up / with ejection pin


	Bite Width	Length of Shaft	
		7" / 18 cm	8" / 21 cm
40° Up-Angle 			
	2 mm	404-540EJ	404-402EJ
	3 mm	404-541EJ	404-404EJ
	4 mm	404-542EJ	404-406EJ
	5 mm	404-543EJ	404-408EJ



Kerrison Punch (Ocular-plasti) 90° Rongeur

4" (10 cm) shaft

16-2100	1.0 mm wide	up
16-2101	1.5 mm wide	up
16-2102	2.0 mm wide	up
16-2106	2.5 mm wide	up
16-2103	3.0 mm wide	up
16-2104	3.0 mm wide	down
16-2105	4.0 mm wide	down
16-2107	4.0 mm wide	up



16-2101



16-2102



16-2103



16-2104



16-2105



16-2106








16-2107

Kerrison Punch (Ocular-plasti) 40° Rongeur

4" (10 cm) shaft

- 16-2096 1.0 mm wide
- 16-2097 1.5 mm wide
- 16-2098 2.0 mm wide
- 16-2099 3.0 mm wide


		Bite Width	
40° Up-Angle		1 mm	16-2096
		1.5 mm	16-2097
		2 mm	16-2098
		3 mm	16-2099



Sontec Kerrison Rongeur

40° / easy clean / 8" (20.3 cm)

- 404-400E 1 mm
- 404-402E 2 mm with ejector pin
- 404-404E 3 mm with ejector pin
- 404-406E 4 mm with ejector pin
- 404-408E 5 mm with ejector pin

	Bite Width	With Pin	Without Pin
40° Up-Angle 	1 mm	NA	404-400E
	2 mm	404-402E	NA
	3 mm	404-404E	NA
	4 mm	404-406E	NA
	5 mm	404-408E	NA



404-402E

404-400E

Rongeurs Tray

Anodized aluminum tray / 9" x 13" x 6"

18-3562



18-3562



Care and Maintenance for the Universal Rotating 360° Kerrison Rongeur Shafts & Handles

The Universal Rotating 360° Kerrison Rongeur


shafts are comprised of four parts:

A: outside shaft with frontal jaw blade

B: Inner shaft with threading at the tip

C: Back side blade/ jaw with screw port at the end


D: spring



Slide inner shaft (B) through outer shaft (A).
Then screw jaw (C) into inner shaft (B).




Screw inner shaft (B) and jaw (C) completely together (only
after Inner shaft (B) and outer shaft (A) are together.)




To put the jaws together, fit top jaw into the grooves
of the outside shaft (A) and slide into place.

Use the grooves on the outside shaft (A) and the back
side jaw/ blade (C) to fit the two parts into position.



 This is the point where the jaw/ blade (C) and inner shaft (B)
are screwed together. You must make sure the two parts are
fully screwed together for the Rongeur to work properly.


INSERTING SHAFT



IMPORTANT:

WHEN INSERTING NEW SHAFT INTO NEW HANDLE, YOU HAVE TO MAKE SURE THAT HANDLE IS IN **CLOSED** POSITION AND MOVING HANDLE FREE TO MOVE. OTHERWISE, THE NEW FORCE-OPEN MECHANISM WILL NOT FUNCTION!

RELEASING SHAFT



Simultaneously with pressing the button pull the lower shaft forward to detach it from the handle.

Click Kerrison punches are suitable for all kinds of shafts: Straight, Bayonet, Curved up and down. You do not need to keep a large stock of Kerrisons, only the shafts are necessary.

This new kind of Kerrison's approach, besides significant financial savings simplifies the entire working process. Your productivity will not be affected by a very short training period. The costs of migration to this new product are very

low and the positive results are almost immediately visible.

For detailed information regarding training and available additional instruction see the contacts on the back cover of this brochure.



Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.
7248 South Tucson Way
Centennial, Colorado 80112 USA
Toll-free / 800.821.7496
Local / 303.790.9411
Fax / 303.792.2606

SontecInstruments.com