
 Sontec Instruments®

SPINE LINE

M.I.S. Straight Tube & Split Tube Sets


## Spine Line MIS Basic Straight Tube Set

SI-TS1	set		
SI-CASETI	case TI one level case		
SI-D108	dilator #1	8 mm	silver
SI-D213	dilator #2	13 mm	green
SI-D318	dilator #3	18 mm	blue
SI-D422	dilator #4	22 mm	red
SI-D526	dilator #5	26 mm	violet
SI-1850T	tube retractor	18 mm x 50 mm	
SI-1860T	tube retractor	18 mm x 60 mm	
SI-1870T	tube retractor	18 mm x 70 mm	
SI-1880T	tube retractor	18 mm x 80 mm	
SI-2250T	tube retractor	22 mm x 50 mm	
SI-2260T	tube retractor	22 mm x 60 mm	
SI-2270T	tube retractor	22 mm x 70 mm	
SI-2280T	tube retractor	22 mm x 80 mm	
SI-2650T	tube retractor	26 mm x 50 mm	
SI-2660T	tube retractor	26 mm x 60 mm	
SI-2670T	tube retractor	26 mm x 70 mm	
SI-2680T	tube retractor	26 mm x 80 mm	



## Spine Line MIS Extended Straight Tube Set

SI-TS2	set		
SI-CASETII	case TII two level case		
SI-D108	dilator #1	8 mm	silver
SI-D213	dilator #2	13 mm	green
SI-D318	dilator #3	18 mm	blue
SI-D422	dilator #4	22 mm	red
SI-D526	dilator #5	26 mm	violet
SI-1840T	tube retractor	18 mm x 40 mm	
SI-1850T	tube retractor	18 mm x 50 mm	
SI-1860T	tube retractor	18 mm x 60 mm	
SI-1870T	tube retractor	18 mm x 70 mm	
SI-1880T	tube retractor	18 mm x 80 mm	
SI-1890T	tube retractor	18 mm x 90 mm	
SI-18100T	tube retractor	18 mm x 100 mm	
SI-18110T	tube retractor	18 mm x 110 mm	
SI-18120T	tube retractor	18 mm x 120 mm	
SI-2240T	tube retractor	22 mm x 40 mm	
SI-2290T	tube retractor	22 mm x 90 mm	
SI-22100T	tube retractor	22 mm x 100 mm	
SI-22110T	tube retractor	22 mm x 110 mm	
SI-22120T	tube retractor	22 mm x 120 mm	
SI-2640T	tube retractor	26 mm x 40 mm	
SI-2650T	tube retractor	26 mm x 50 mm	
SI-2660T	tube retractor	26 mm x 60 mm	
SI-2670T	tube retractor	26 mm x 70 mm	
SI-2680T	tube retractor	26 mm x 80 mm	
SI-2690T	tube retractor	26 mm x 90 mm	
SI-26100T	tube retractor	26 mm x 100 mm	
SI-26110T	tube retractor	22 mm x 110 mm	
SI-26120T	tube retractor	26 mm x 120 mm	



beveled tubes  
in all sizes



22 mm diameter



18 mm diameter




26 mm diameter




## Spine Line MIS Basic Split Tube Set

SI-ST1	set
SI-CASETI	case TI one level case
SI-D108	dilator #1 8 mm silver
SI-D213	dilator #2 13 mm green
SI-D318	dilator #3 18 mm blue
SI-D422	dilator #4 22 mm red
SI-1850MB	18 mm x 50 mm retractor blades (R/L)
SI-1860MB	18 mm x 60 mm retractor blades (R/L)
SI-1870MB	18 mm x 70 mm retractor blades (R/L)
SI-1880MB	18 mm x 80 mm retractor blades (R/L)
SI-2250MB	22 mm x 50 mm retractor blades (R/L)
SI-2260MB	22 mm x 60 mm retractor blades (R/L)
SI-2270MB	22 mm x 70 mm retractor blades (R/L)
SI-2280MB	22 mm x 80 mm retractor blades (R/L)
SI-0125	top-loading floating frame 240 mm rack with pivoting arms




## Spine Line MIS Extended Split Tube Set

SI-ST2	set
SI-CASETII	case TII two level case
SI-7190S	hook straight 90 mm
SI-71100S	hook straight 100 mm
SI-71120S	hook straight 120 mm
SI-1440SS	micro teeth blade 40 mm
SI-1450SS	micro teeth blade 50 mm
SI-2365AS	angled teeth blade 65 mm
SI-2375AS	angled teeth blade 75 mm
SI-2385AS	angled teeth blade 85 mm
SI-23100AS	angled teeth blade 100 mm
SI-23120AS	angled teeth blade 120 mm




18 & 22 mm diameters




new pivoting frame: up to 40 degrees of distal blade rotation

## Optional Items

---



**Light Cable Ports**  
Offer superior directional light



**Light Cable Ports**  
Offer superior directional light



**Articulating Arm**



**Light Source**



**Sontec Instruments®**

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.  
7248 South Tucson Way  
Centennial, Colorado 80112 USA  
Toll-free / 800.821.7496  
Local / 303.790.9411  
Fax / 303.792.2606

[SontecInstruments.com](http://SontecInstruments.com)