




 Sontec Instruments®

Micro V.A.T.S. Uniportal System

Micro V.A.T.S. Uniportal System

2.7 mm


- | | | |
|--------------|---|--|
| 20903-17 | 1 | MICTEC flexible trocar set Ø 3.5 mm, black plastic fixation screw, 70 mm working length |
| 20903-17-310 | 2 | MICTEC flexible trocar set Ø 3.5 mm, black plastic fixation screw, with sealing cap, 70 mm working length |
| 27913-35 | 1 | MICTEC trocar set pyramidal tip with manual folding valve, with insufflation tap, Ø 3.5 mm, 95 mm working length, with smooth sleeve |



Micro V.A.T.S. Uniportal System

2.7 mm

27971-45	1	MICTEC fenestrated forceps, 30° curved, down, Y-grip with clamp ratchet, with cleaning flush port, Ø 2,7 mm, 300 mm working length
27971-90	1	MICTEC atraumatic dressing forceps, 30° curved down, Y-grip without ratchet, with cleaning flush port, Ø 2.7 mm, 300 mm working length
27973-36	1	MICTEC dissection clamp, 30° curved down, Y-grip without ratchet, with cleaning flush port, Ø 2.7 mm, 300 mm working length
27975-37	1	MICTEC scissors curved, 30° curved down, Y-grip, with cleaning flush port, Ø 2.7 mm, 300 mm working length
27976-36	1	MICTEC biopsy forceps, 30° curved down, Y-grip without ratchet, with cleaning flush port, Ø 2,7 mm, 300 mm working length



Micro V.A.T.S. Uniportal System

2.7 mm


26943-21	1	Suction-Irrigation handle with 3-way-tap large adapter Luer-Lock-male
27983-30	1	MICTEC suction and irrigation tube, Ø 2.7 mm, 30° curved down, 300 mm working length
27948-30	1	MICTEC hook electrode, insulated, Ø 2.7 mm, 30° curved up, 300 mm working length
27948-32	1	MICTEC flat electrode, insulated, Ø 3.2 mm, 30° curved up, 300 mm working length
26944-03	1	MICTEC HF-cable with hexagon inside 4m long; connection: Valleylab



Micro V.A.T.S. Uniportal System

2.7 mm


70373-58	1	Perforated Sterilizing Tray with Lid and Silicone Holders 545 x 255 x 95 mm
70373-58-062	1	Perforated Sterilizing Inlay for Tray 360 x 230 x 33.5 mm




70373-58-062



70373-58





Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.
7248 South Tucson Way
Centennial, Colorado 80112 USA
Toll-free / 800.821.7496
Local / 303.790.9411
Fax / 303.792.2606

SontecInstruments.com